

# **CRYO-FLUORESCENCE TOMOGRAPHY**

### Transformative 3D Imaging to Monitor Drug PK/PD

Cryo-Fluorescence Tomography (CFT) is a transformative 3D approach to imaging drug distribution and protein expression in whole-animals. EMIT Imaging offers both instrumentation and services via our platform, Xerra<sup>™</sup>, a high resolution and high sensitivity preclinical imaging tool designed to advance discoveries in biological and drug research. See what you're missing with standard 2D *in vivo* fluorescence imaging techniques and utilize CFT to better visualize:

- Whole-body drug distribution and delivery
- Screening candidate drugs and delivery systems
- Whole-body therapeutic protein expression
- Multiplexed co-localization of drug with targets
- On and off target tissue identification



# **HOW CFT WORKS**

## 1.PREP

Cryo-preserve & embed the sample

## 2. IMAGE

High resolution, high sensitivity white light & flourescence imaging

## **3. SECTION**

Remove 20–55  $\mu m$  from the block surface

## OUTPUT

3D image stacks









HIGH RESOLUTION Resolution down to 20 µm HIGH SENSITIVITY nM sensitivity, comparable to nuclear medicine **MULTIPLEXING** 6 lasers and 7 filters for multiplexed applications

AUTOMATION Xerra is a fully automated



## **APPLICATIONS**

### Fibroblast-activation protein (FAP): Whole-Body Drug Biodistribution



3D Fluorescence Image Anatomical White-light Fluorescence Image

Fluorescence

Anatomical +

Copyright 2023 Ratio Therapeutics Inc., used with permission.

FINDINGS: Data confirms drug preferential uptake in the tumor with minimal uptake in other organs/tissues

#### **DETAILS:**

- 1µg dose ofZ W800-1-labeled FAP-targeting small molecule
- U87 xenograft
- 4hrs post iv administration
- Imaged at 35 μm

### **AAV-mediated Protein Expression Following IV Administration**



AAV9 Control



**FINDINGS:** CFT uniquely achieves the required combination of resolution and sensitivity in a whole animal

#### **DETAILS:**

- GFP expression at day 15 following IV injection
- No expression in the control
- Liver is major organ where GFP expression occurs
- Note signal in dorsal root ganglion (DRG) and heart

AAV9-GFP



## **APPLICATIONS**

### Whole-Body Distribution of a Labeled ASO Following IT Injection



2D In Vivo Fluorescence



3D CFT

**FINDINGS:** CFT is superior to surface fluorescence due to its depth limitations (signal loss) and poor resolution

#### **DETAILS:**

- Cy5-labeled ASO at day 15 following IT injection
- Full visualization of spinal column into the cisterna magna
- CFT detects bioD and clearance: GI, lympathic, nasal turbinate, etc.
- CFT is 3D and high resolution:
  <50 μm versus >1mm in optical and nuclear medicine

#### **Cell-specific Fluorescence Reporter Imaging**





#### **DETAILS:**

GAD67-GFP (line G42) mice selectively express enhanced green fluorescent protein (EGFP) in the parvalbumin (Pv)-expressing subclass of basket interneurons (soma, dendrites, and axons) and also in putative presynaptic boutons.



## **XERRA<sup>™</sup> CFT IMAGING SYSTEM**

#### **INSTALL XERRA IN YOUR LAB**

- Xerra automates the CFT workflow
- Anatomical white-light images co-registered with fluorescence images, capable of multiplexing
- 5 magnifications: 20-55 µm pixel resolution
- 6 excitation lasers: 470 to 780 nm
- 7 emission filters: 500 to 850 nm
- Xerra is CE marked

#### **RUN CFT AS A SERVICE**

- EMIT offers services, including imaging and analysis
- Animals/tissues are prepared at the sponsor or CRO
- Shipped frozen to EMIT for blocking and imaging
- Images are reviewed with the sponsor and additional analysis support available



## **FLUOROPHORE CHARACTERIZATION**

Xerra is equipped with laser and filter combinations that match most applications in biological research



The Limit of Detection (LOD) and Limit of Quantification (LOQ) have been evaluated for the following fluorophores:

Fluorophore	LOD (nM)	LOQ (nM)
СуЗ	1.7	5.1
Cy5	2.6	7.7
Cy5.5	4.1	12.3
AF488	4.3	12.9
AF546	6.8	20.5
AF647	0.9	2.9
AF700	3.8	11.7
AF750	3.7	11.3
IRDye 800CW	9.5	28.9
ZW800-1	11.7	35.3

Fluorophore Queue:

- GFP
- mCherry
- tdTomato
- mKate
- FITC
- Texas Red
- AF594

Additional can be added upon request